





## xCELLigence RTCA DP Instrument

*Flexible Real-Time Cell Monitoring*



**For life science research only.  
Not for use in diagnostic procedures.**

# The xCELLigence RTCA DP Instrument

## Flexible Real-Time Cell Monitoring


The RTCA DP Instrument expands the throughput and application options of the xCELLigence Real-Time Cell Analyzer (RTCA) portfolio. Featuring a dual-plate (DP) format, the instrument measures impedance-based signals in both cellular and cell invasion/migration (CIM) assays – without the use of exogenous labels. With outstanding application flexibility, the RTCA DP Instrument supports multiple users performing short-term and long-term experiments.

### Explore the wide range of applications


- Cell invasion and migration assays
- Compound- and cell-mediated cytotoxicity
- Cell adhesion and cell spreading
- Cell proliferation and cell differentiation
- Receptor-mediated signaling
- Virus-mediated cytopathogenicity
- Continuous quality control of cells

The xCELLigence System continuously and non-invasively detects cell responses throughout an experiment, without the use of exogenous labels that can disrupt the natural cell environment.

- **Obtain complete, continuous data profiles** from cell responses generated during *in vitro* experiments (Figure 1).
- **Take advantage of real-time data** to identify optimal time points for downstream assays.
- **Combine real-time monitoring of cellular responses with complementary functional endpoint assays**, and maximize data quality before, during, and after your experiment.



**Figure 1: Reveal cytotoxic effects through continuous monitoring.** HT1080 cells were seeded in an E-Plate at two different densities (5,000 and 10,000 cells) and treated 24 hours later with 12.5 nM Paclitaxel, or DMSO as a control. As shown by the Cell Index profile, which reflects cell adherence, the antimitotic effect of Paclitaxel was observed in HT1080 cells that were proliferating, whereas confluent cells showed no response.



### Compact. Convenient. Versatile.

The RTCA DP Instrument consists of two components: the RTCA Control Unit and the RTCA DP Analyzer with three integrated stations for measuring cell responses in parallel or independently.

- Choose from three types of impedance-based 16-well plates:
  - E-Plate 16 and E-Plate VIEW 16 for cellular assays
  - CIM-PLATE 16 for cell invasion/migration assays
- Use all three different plate types in any combination.
- Easily achieve optimal cell culture conditions by placing the RTCA DP Analyzer and plates into standard CO<sub>2</sub> incubators.

## E-Plates for the RTCA DP Instrument

*More Flexibility. More Data. More Insight.*


Obtain detailed information about your cells with the versatile RTCA DP Instrument, which supports up to three plates of any type – E-Plate 16, E-Plate VIEW 16, or CIM-Plate 16 – in any combination. For example, cell invasion/migration assays and cytotoxicity assays or short- and long-term assays may be run simultaneously.

### E-Plate 16 and E-Plate VIEW 16: Cellular Assays in a 16-Well Format



- Quantitatively monitor changes in cell number, cell adhesion, cell viability, and cell morphology.
- Easily add compounds during an experiment.
- Assess short- and long-term cellular effects.
- With the E-Plate VIEW 16, observe measured changes using microscopes.



E-Plate 16



**Figure 2: Easily visualize cells while measuring cell response with xCELLigence System E-Plate VIEW technology.** A modified version of the standard E-Plate 16, the E-Plate VIEW 16 enables image acquisition using microscopes or automated cell-imaging systems. For the modification, four rows of microelectrode sensors were removed in each well to create a window for visualizing cells. Approximately 70% of each well bottom is covered by the microelectrodes, providing cell impedance measurements nearly identical to those obtained with the standard E-Plate 16. Both plate types can be used in parallel.




**Figure 3: Continuously monitor cells and determine optimal time points for assessing cytotoxicity.** Cell proliferation and cell death were continuously monitored using the xCELLigence RTCA DP Instrument. The optimal time points for visual inspection of HeLa cells were determined and images taken 4 and 22 hours after compound treatment using a Z16 Apo Microscope with light base (Leica Microsystems).

## CIM-Plate 16: Quantitative Cell Invasion/Migration Analysis


- Monitor cell invasion and migration continuously in real time over the entire time course of an experiment.
- Eliminate time-consuming manual detection (Figure 4).
- Perform CIM analysis in a convenient one-well system (Figure 5).



CIM-Plate 16



**Figure 4: Quantitatively measure the rate and onset of invasion while concurrently assessing migration.** HT1080 cells ( $2 \times 10^4$ ) were seeded in the upper chamber of CIM-Plate wells coated with varying dilutions of Matrigel, or in wells with no coating. Serum was added to the lower chamber of selected wells as a chemoattractant. Invasion was observed and migration monitored continuously over a 70-hour period. All serum-starved samples resulted in base-line Cell Index levels, indicating the absence of invasion/migration, while those wells with chemoattractant induced migration.



**Figure 5: Analyze invasion/migration in real time with the CIM-Plate 16.** The plate features two separable sections for ease of experimental setup. Cells seeded in the upper chamber move through the microporous membrane into the lower chamber that contains a chemoattractant. Cells adhering to the microelectrode sensors lead to an increase in impedance, which is measured in real time by the RTCA DP Instrument.

## E-Plate Insert 16: Co-Culture in Real-Time


- Continuously monitor indirect cell-cell interactions.
- Assess short- and long-term cell response without labor-intensive labeling and microscopy.
- Co-culture different cell types under physiological conditions for a broad range of applications, including:

**Cancer Research:** Assess paracrine stimulation of cancer cell proliferation by fibroblasts.

**Immunology:** Investigate immune cell interactions.


**Stem Cell Research:** Monitor proliferation and differentiation in the presence of stimulation cells.

**Toxicology:** Determine cytotoxicity of agents and assess effects of cytokine release.



**Figure 6. Real-time monitoring of co-culture-induced proliferation stimulation and its inhibition using the E-Plate Insert.** Intercellular interactions play an important role in normal cell development and tumorigenesis. Results show that the proliferation of hormone-responsive tumor cells is likely mediated by hormones and growth factors exchanged between the two cell populations separated by the E-Plate Insert.

Elevated T47D cell proliferation on the E-Plate (green trace ■) was induced by hormone secretion of H295R cells in the insert, and inhibited by the hormone synthesis inhibitor Prochloraz (blue trace ■). Incubation of T47D cells with only the E-Plate Insert did not affect proliferation (red trace ■).



# Selected Publications for the RTCA DP Instrument

## 1. Cell Invasion and Migration

### **MicroRNA-200c Represses Migration and Invasion of Breast Cancer Cells by Targeting Actin-Regulatory Proteins FHOD1 and PPM1Ferences.**

Jurmeister S, Baumann M, Balwierz A, Keklikoglou I, Ward A, Uhlmann S, Zhang JD, Wiemann S, Sahin O. *Mol Cell Biol.* 2012; 32(3):633–651.

### **c-Myb regulates matrix metalloproteinases 1/9, and cathepsin D: implications for matrix-dependent breast cancer cell invasion and metastasis.**

Knopfová L, Beneš P, Pekarčíková L, Hermanová M, Masařík M, Pernicová Z, Souček K, Smarda J. *Mol Cancer.* 2012; 11:15.

### **Comparative Analysis of Dynamic Cell Viability, Migration and Invasion Assessments by Novel Real-Time Technology and Classic Endpoint Assays.**

Limame R, Wouters A, Pauwels B, Fransen E, Peeters M, Lardon F, De Wever O, Pauwels P. *PLoS One.* 2012; 7(10): e46536.

## 2. Compound-mediated Cytotoxicity/Apoptosis

### **Screening and identification of small molecule compounds perturbing mitosis using time-dependent cellular response profiles.**

Ke N, Xi B, Ye P, Xu W, Zheng M, Mao L, Wu MJ, Zhu J, Wu J, Zhang W, Zhang J, Irelan J, Wang X, Xu X, Abassi YA. *Anal Chem.* 2010; 82(15):6495-503.

### **Kinetic cell-based morphological screening: prediction of mechanism of compound action and off-target effects.**

Abassi YA, Xi B, Zhang W, Ye P, Kirstein SL, Gaylord MR, Feinstein SC, Wang X, Xu X. *Chem Biol.* 2009; 16(7):712-23.

## 3. Cell-mediated Cytotoxicity

### **Real-time profiling of NK cell killing of human astrocytes using xCELLigence technology.**

Moodley K, Angel CE, Glass M, Graham ES. *J Neurosci Methods.* 2011; 200(2): 173-180.

### **Unique functional status of natural killer cells in metastatic stage IV melanoma patients and its modulation by chemotherapy.**

Fregni G, Perier A, Pittari G, Jacobelli S, Sastre X, Gervois N, Allard M, Bercovici N, Avril MF, Caignard A. *Clin Cancer Res.* 2011; 17(9): 2628–37.

## 4. Cell Adhesion and Cell Spreading

### **A role for adhesion and degranulation-promoting adapter protein in collagen-induced platelet activation mediated via integrin $\alpha 2 \beta 1$ .**

Jarvis GE, Bihan D, Hamaia S, Pugh N, Ghevaert CJ, Pearce AC, Hughes CE, Watson SP, Ware J, Rudd CE, Farndale RW. *Journal of Thromb Haemost.* 2012; 10(2): 268–277.

### **Dynamic monitoring of cell adhesion and spreading on microelectronic sensor arrays.**

Atienza JM, Zhu J, Wang X, Xu X, Abassi Y. *J Biomol Screen.* 2005; 10(8): 795-805.

## Selected Publications continued

### 5. Receptor-mediated Signaling

**Impedance responses reveal  $\beta_2$ -adrenergic receptor signaling pluridimensionality and allow classification of ligands with distinct signaling profiles.**

Stallaert W, Dorn JF, van der Westhuizen E, Audet M, Bouvier M.  
*PLoS One*. 2012; 7(1): e29420.

**Label-free impedance responses of endogenous and synthetic chemokine receptor CXCR3 agonists correlate with Gi-protein pathway activation.**

Watts AO, Scholten DJ, Heitman LH, Vischer HF, Leurs R.  
*Biochem Biophys Res Commun*. 2012; 419(2):412-8.

**Impedance measurement: A new method to detect ligand-biased receptor signaling.**

Kammermann M, Denelavas A, Imbach A, Grether U, Dehmlow H, Apfel CM, Hertel C.  
*Biochem Biophys Res Commun*. 2011; 412(3): 419-424.

### 6. Virus-mediated Cytopathogenicity

**Novel, real-time cell analysis for measuring viral cytopathogenesis and the efficacy of neutralizing antibodies to the 2009 influenza A (H1N1) virus.**

Tian D, Zhang W, He J, Liu Y, Song Z, Zhou Z, Zheng M, Hu Y.  
*PloS One*. 2012; 7(2):e31965.

**Real-time monitoring of flavivirus induced cytopathogenesis using cell electric impedance technology.**

Fang Y, Ye P, Wang X, Xu X, Reisen W.  
*J Virol Methods*. 2011; 173(2):251-8.

### 7. Quality of Control of Cells

**Rapid and quantitative assessment of cell quality, identity, and functionality for cell-based assays using real-time cellular analysis.**

Irelan JT, Wu MJ, Morgan J, Ke N, Xi B, Wang X, Xu X, Abassi YA.  
*J Biomol Screen*. 2011; 16(3):313-22.

**Live cell quality control and utility of real-time cell electronic sensing for assay development.**

Kirstein SL, Atienza JM, Xi B, Zhu J, Yu N, Wang X, Xu X, Abassi YA.  
*Assay Drug Dev Technol*. 2006; 4(5):545-53.


### 8. Endothelial Barrier Function

**An inverted blood-brain barrier model that permits interactions between glia and inflammatory stimuli.**

Sansing HA, Renner NA, MacLean AG.  
*J Neurosci Methods*. 2012; 207(1):91-6.

**A dynamic real-time method for monitoring epithelial barrier function in vitro.**

Sun M, Fu H, Cheng H, Cao Q, Zhao Y, Mou X, Zhang X, Liu X, Ke Y.  
*Anal Biochem*. 2012; 425(2):96-103.



## Ordering Information for xCELLigence RTCA DP System

Product	Cat. No.	Pack Size
<b>xCELLigence RTCA DP Instrument</b>	<b>00380601050</b>	<b>1 Bundled Package</b>
RTCA DP Analyzer	05469759001	1 Instrument
RTCA Control Unit	05454417001	1 Notebook PC
<b>E-Plate 16</b>	05469830001	6 Plates
	05469813001	6 x 6 Plates
<b>E-Plate VIEW 16</b>	06324738001	6 Plates
	06324746001	6 x 6 Plates
<b>E-Plate Insert 16</b>	06465382001	1 x 6 Devices (6 16-Well Inserts)
<b>CIM-Plate 16</b>	05665817001	6 Plates
	05665825001	6 x 6 Plates
<b>CIM-Plate 16, Assembly Tool</b>	05665841001	1 Assembly Tool

Learn more about the enabling technology of the xCELLigence System and its broad range of applications at [www.aceabio.com](http://www.aceabio.com)

**For life science research only.  
Not for use in diagnostic procedures.**

XCELLIGENCE, E-PLATE, CIM-PLATE, and ACEA BIOSCIENCES are registered trademarks of ACEA Biosciences, Inc. in the US and other countries.

All other product names and trademarks are the property of their respective owners.

### Published by

ACEA Biosciences, Inc.  
6779 Mesa Ridge Road Ste. 100  
San Diego, CA 92121  
U.S.A.

[www.aceabio.com](http://www.aceabio.com)

© 2013 ACEA Biosciences, Inc.  
All rights reserved.